





Connecting OMRON Dual Channel Series 7 Model BP7350 to OMRON App and syncing with COACH

Downloading the OMRON Connect app and BP Cuff pairing

1. Download OMRON Connect Mobile App on App store or Google Play.




2. Make an OMRON account.
3. To connect your Cuff to your Device, make sure bluetooth is turned ON. Open the App and click on .



4. Click on “profile” at the bottom of screen




5. Click on “Connected devices”



6. Click on “Add new device”



7. On the next screen select “All OMRON Blood Pressure devices” and tap on “Next: Pairing” at the bottom of screen.
8. To put your BP monitor in pairing mode press the bluetooth button for 3-5 seconds. On your Series 7 Model BP7350 it will be located on back side as shown below. Once in pairing mode tap the “Next: Begin pairing” button at bottom of your app.



9. Select the device you want to pair. Once paired your monitor screen will display "OK".



10. Now you are ready to transfer your blood pressure readings to the OMRON Connect app from your BP Cuff.



Syncing OMRON Connect to COACH App

1. On the COACH home page click on “Click here”

COACH Home Goals Home BP Readings My Medications Resources Refresh Log Out

Andrew Fhir - 42 years Male Clear Supplemental Data

High Blood Pressure Control

Your blood pressure:

STOP

Most Recent BP:
98 systolic
55 diastolic

Warning: Your BP is very low. Take action below.

Click here to authenticate and synchronize with Omron.

Your Current Blood Pressure Goal: Below 140 / 90 (update)

Recommendations Monitoring

2. You will be directed to OMRON login page. Use your email and password for OMRON Connect app.



Login

Email Address

Password

Login

[Forgot password?](#)



3. Check the boxes for date you wish to share with COACH and click “Allow”

OMRON

COACH
is requesting your permission

Check the permissions you wish to grant.

Personal Information

☒ **Your e-mail**(required)

Personal Information

Access your blood pressure readings
☒ This application is requesting access to your blood pressure related data stored by Omron

Access your fitness metrics
☒ This application is requesting access to your fitness related data stored by Omron

Allow

Do Not Allow



- Once you're back in COACH, take a quick look at the home page. You'll see the date and time your data was last synced. This ensures all your information is up-to-date.

HomeGoalsHome BP ReadingsMy MedicationsResources

RefreshLog Out

Andrew Fhir - 42 years Male

Clear Supplemental Data

High Blood Pressure Control

Your blood pressure:

Most Recent BP:

98 systolic

55 diastolic

Warning: Your BP is very low. Take action below.

Omron data last synchronized Friday, December 15, 2023 at 11:59 AM

Your Current Blood Pressure Goal: **Below 140 / 90** [\(update\)](#)

Recommendations

Monitoring